

# Resus Checklist - Bay 2

Week commencing:

	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
Initials of person checked this bay							
<b>Left Side of Bay</b>							
Adult bag Valve Mask (drip stand)							
Oxygen ports intact and working							
Suction unit working with tubing and Yankeur (in packaging) connected							
NRB mask on top of Intubation box next to O2 port							
X2 O2 cylinders (FULL) underneath Trolley							
X1 Spare Suction Canister (on 2 <sup>nd</sup> Shelve)							
Area cleaned (Green sticker applied)							
<b>Left Side Top Draw</b>							
Nasopharyngeal x2 Sizes - 6.0 + 7.0							
Oropharyngeal x2 Size - 2, 3, 4							
Vent Masks Sizes - 3, 4, 5							
LMA / IGEL Sizes 3, 4, 5							
Aqua Gel x3							

<b>Left Side Bottom Draw</b>							
Oxylog 3000 ventilation tubing							
Adult Breathing Circuit							
<b>Right Side Of Bay</b>							
Observations monitor Cleaned and working with all leads present. (ECG, BP, Sats) c							
Defibrillator Machine plugged in and daily op checked with ECG leads present							
Adult Defibrillator pads X2							
ECG Electrodes x2							
Stethoscope (Right Drip Stand)							
Graesby Pump (Plugged In)							
Aguilla Pump (Plugged In)							
500 ml / 1000ml Fluid Pressure Bags (Right Drip Stand)							
Area Cleaned (Green Sticker applied)							
<b>Right Side Draw</b>							
Blood Pressure Cuff (Small, Medium, Large)							
Rectal Temperature Sensor							
Gauze x2							
Small Bandages X2							
Atropine Minijet x2							
Adrenaline Minijet x6							
Amiodarone Minijet x2							
Spare Defib Paper x1							

Entonox Mouth Pieces x2							
Razors x2							
<b>Back Of Bay</b>							
Suction Catheters x2 Sizes 6 – 16							
Yankeur's x4							
Bougie Sizes 5G, 10G, 15G							
Stylet Sizes 2.0, 4.0							
Extra long Stylet							
Ventilator (plugged in) with Oxylog Tubing							
X2 Blue Emergency Drug Boxes (sealed and in date)							
Trolley cleaned with Tristel (Green sticker applied)							
<b>Intubation box</b>							
ET Tubes cuffed x2 sizes 5.0 – 9.0							
Laryngoscope Handle x1 (check working)							
Laryngoscope blades – MAC 2,3,4 + LT3							
Adult Magills forceps							
10ml syringe							
Sterile Scissors							
Aqua gel x3							
ET tube tie / Ribbon							
Philips Filter line							
Catheter mount (filter and T piece attached)							

# Bay 2 Additional checks

	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
Initial of person checking							
Oxygen and Hygiene Trolley stocked Cleaned (Green sticker applied)							
Arterial Line Box Checked / Sealed							
CVP Box Checked / Sealed							
Chest Drain Box Checked / Sealed							
Bladder Scanner cleaned and stocked with Aqua Gel (Green Sticker applied)							
NIV Trolley stocked and cleaned (Green Sticker applied)							
Catheter trolley stocked and cleaned (Green sticker applied)							