A&E Daily Cubicle check list

	Oxygen cylinder under trolley - Working & trolley cleaned with Tristell	Wall oxygen - working and NRB bag with port	Suction – working, clean and filter in date.	Nurse call buzzer – Check its working, not damaged.	Obs machine, leads and BP cuffs cleaned – Damp clean with Tristell	Opthalmoscope/ Otoscope working Damp clean with Tristell	Shelf wiped down & Restock vomit bowls	Boards – Up to date.
Cubicle 1								
Cubicle 2								
Cubicle 3								
Cubicle 4								
Cubicle 5								

Remember to green sticker all areas that have been cleaned

Date:

Person checking:

ESCALATE TO NURSE IN CHARGE ANY ISSUES & MAKE SURE IT LOGGED AS A JOB.