A&E Daily Cubicle check list

	Oxygen cylinder under trolley - Working & trolley cleaned with Tristell	Wall oxygen - working and NRB bag with port	Suction – working, clean and filter in date.	Nurse call buzzer – Check its working, not damaged.	Obs machine, leads and BP cuffs cleaned – Damp clean with Tristell	Opthalmoscope/ Otoscope working	Shelf wiped down & Restock vomit bowls	Boards – Up to date.
Cubicle 6								
Cubicle 7								
Cubicle 8								
Cubicle 9								
Cubicle 10								
Cubicle Paeds								

Remember to green sticker all areas that have been cleaned

D	2	ŧ	Δ	-
$\boldsymbol{\mathcal{L}}$	u		v	•

Person checking:

ESCALATE TO NURSE IN CHARGE ANY ISSUES & IF REQUIRED MAKE SURE IT IS LOGGED AS A JOB.