

**Resus Paeds Bay Checklist**

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thurs</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>Left Side</b>							
FACUDA <b>**Plugged in and Working**</b> , SpO2, ECG BP & Capnography leads							
Philips defib <b>**Plugged in &amp; Daily Op Check**</b> , Test block & ECG leads							
Adult Defib Pads							
Paeds Defib Pads							
ECG Electrodes x1 pack							
X2 Red drug boxes (in date & seal intact)							
<b>Right Side</b>							
Disposable Kidney Dishes <b>**No more than 10**</b>							
Capnography lead							
BP Cuffs: Adult, Small Adult, Child, Small child, Infant, Small Infant & Disposable adult, X1 each <b>**Cleaned daily with Tristel**</b>							
ET Tube Tie (approx 1 metre)							
O2 point working with Paed & Adult NRB mask							
Sharps Bin (Replace if needed)							
Pressure bags: 1000ml, 500ml 1 of each							

Paediatric Resus Checklist

Otoscope (with covers) & Ophthalmoscope **Plugged in both working**							
Stethoscope (Paediatric)							
Breathing System for BabyPac Vent Ref 6203							
PneuPac BabyPac 100 Ventilato							
Adult Mapelson C Circuit ***Un-opened still in Packaging*** Ref 2102000							
Child Mapelson F Circuit ***Un-opened still in Packaging*** Ref 2121000							
Breathing System Neonatal Ref 6362000							
Suction unit (Working) with tubing (still in packaging) ** Does filter need changing ???**							
TransWarmer x1							
Adult & Paeds Yankauer x3							
Blue drug box (in date & seal intact)							
Trolley Cleaned with Tristel Including underneath.							
Blood Trolleys Cleaned with Tristel							
Nurse station clean & clutter Free							

Paediatric Resus Checklist

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thur</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>Paediatric Circulation Trolley</b>							
<b>IV Access Drawer 1</b>							
Hypoglycaemia Pack							
IV Packs x3							
Jelco x5 Yellow & x3 Blue							
Steristrips x4 Small							
Tegaderm x2							
10 x 10 Gauze x3 packs							
5cm K-Band x2							
Micropore roll x1							
5ml syringes x3							
10ml syringes x2							
Red filter needle x4							
Blue needles x3							
Blood Culture Bottle (Pink Top) x1							
Blood bottles FBC, U&E, Glucose, X- Match ***x5 of Each***							
Blue & Yellow 'Stabbers' x2 each							
Tube of Yellow Parafin							
Capillary tubes, Caps & Clot Catchers							
1 x Each size hand cradles							

Paediatric Resus Checklist

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thurs</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>Paediatric Circulation Trolley</b>							
<b>IV Access Drawer 2</b>							
Hypoglycaemia Pack (separate checklist inside)							
IV Packs x3							
Jelco x5 Yellow & x3 Blue							
Steristrips x4 Small							
Tegaderm x2							
10 x 10 Gauze x3 packs							
5cm K-Band x2							
Micropore roll x1							
5ml syringes x3							
10ml syringes x2							
Red filter needle x4							
Blue needles x3							
Blood Culture Bottle (Pink Top) x1							
Blood bottles FBC, U&E, Glucose, X-Match ***x5 of Each***							
Blue & Yellow 'Stabbers' x2 each							
Tube of Yellow Parafin							
Capillary tubes, Caps & Clot Catchers							
1 x Each size hand cradles							

Paediatric Resus Checklist

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thurs</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>Paediatric Circulation Trolley</b>							
<b>Fluid Drawer</b>							
2 x 100ml 0.9% Normal Saline							
2 x 500ml 0.9% Normal Saline							
2 x 500ml 0.9% Normal Saline & 5% Glucose							
1 x 500ml 10% Glucose							
1 x Human Albumen Solution							
2 x Giving sets for Alaris SE pump							
1 x Blood giving set for Alaris SE pump							
2x Single catheter extension sets							
2 x 3-Way Taps							
6 x Red Drawing up Needles							
1 x 50ml Luer-Lok Syringe							
6 x 20ml Luer-Lok Syringe							
6 x 10ml Luer-Lok Syringe							
Mini Spike for drawing up Albumen							

Paediatric Resus Checklist

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thur</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>Paediatric Circulation Trolley</b>							
<b>Intraosseous Access Drawer</b>							
EZ-IO Intraosseous Drill							
EZ-IO Needle 3-39KG (Pink) x3							
EZ-IO Needle 40 KG or over (Blue) x3							
EZ-Connect Fluid transfer set x6							
Cook Medical Intraosseous Needle 15.5G x3							
Chloraprep Sponge x2							
Syringes 50ml x1, 10ml x2 & 1ml x2							
Red Drawing Up Needle x2							
Orange needle x2							
0.9% NaCl Flush 10ml x2							
EZ-IO Instruction CD							

Paediatric Resus Checklist

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thur</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>Catheterisation Grab Box (Paed)</b>							
Newcastle Urine Collection Pad x2							
2 x Red Top Urine Bottles							
4 x U-Bags							
Urometer							
50ml Catheter Tip Syringe x1							
Catheters Size 6, 8 &10 x2 each							
Dressing Pack x1							
1 x 10ml ampules Sterile Water							
1 x 10ml Syringe							
1 x Red Drawing up Needle							
GlycoBloc x1							
Instilagel x1							
S/M/L Sterile Gloves x1 each							

Paediatric Resus Checklist

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thur</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>Paediatric CVP Drawer</b>							
Sterile Gown & Sterile Drape x1 each							
Sterile Gloves S/M/L x1 each							
Chloraprep 1.5ml applicator x1							
2.5 fr Careflow Peripheral Venous Catheter x1							
Secalon 16G central venous catheter x2							
Portacath Needle 19G 3/4"							
Large Tegaderm x2							
10ml & 5ml Syringe x1 each							
5ml ampules Lidocaine 1% x2							
2.0 Sofsilik Suture x1							
10ml Ampules 0.9% sodium Chloride x2							
Orange needle x1							
Red filter needle x1							
Sterile Disposable Scissors x 1							
Blue needle free valve x1							



Paediatric Resus Checklist

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thur</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>Paediatric Airway Trolley</b>							
<b>O2 Masks &amp; Nebs</b>							
Adult NRB mask x2							
Paed NRB mask x2							
Adult Neb x1							
Paed Neb x2							
Nasal Specs Neonatal x1							
Nasal Specs Paed x1							
Nasal Specs Adult x1							
<b>Face Masks</b>							
Face Mask size 0 (round) x2							
Face Mask size 2 x2							
Face Mask size 3 x2							
Face Mask size 4 x1							
Face Mask size 5 x1							
<b>Oropharyngeal &amp; Nasopharyngeal</b>							
OPA size 000 x2							
OPA size 00 x2							
OPA size 0 x2							
OPA size 1 x2							
OPA size 1.5 x2							
OPA size 2 x4							
OPA size 3 x3							
OPA size 4 x2							
NPA size 3 x1							
NPA size 4 x1							
NPA size 5 x1							
NPA size 6 x1							
NPA size 7 x1							
NPA size 8 x1							
Tongue Depressor x2							

Paediatric Resus Checklist

Day Insert Date	Mon	Tue	Wed	Thur	Fri	Sat	Sun
Staff Member							
<b>Paediatric Airway Trolley</b>							
<b>Orogastric tubes</b>							
OG tube size 6 x1							
OG tube size 8 x1							
OG tube size 10 x1							
OG tube size 12 x1							
OG tube size 14 x1							
OG tube size 18 x1							
3 way taps x2							
Duoderm x1 sheet							
Tegaderm Dressing 5							
Pentorch x1							
Syringe 10ml x1							
Syringe 20ml x1							
Neofit x1							
Bile Bag x1							
<b>Bag Valve &amp; Mask</b>							
Adult BVM x1							
Child BVM x1							
Infant BVM x1							
Mapelson F - Paediatric (Ref: 2121000) x1							
<b>Laryngoscope</b>							
Laryngoscope handle Large x2							
Laryngoscope handle Small x2							
Laryngoscope Blade Mac 2 x2							
Laryngoscope Blade Mac 3 x2							
Laryngoscope Blade Mac 4 x2							
Laryngoscope Blade - Miller size 00 x1							
Laryngoscope Blade - Miller size 0 x1							
Laryngoscope Blade - Miller size 1 x1							
Batteries size C x2							
Batteries size AA x2							

Paediatric Resus Checklist

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thur</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member							
<b>Paediatric Airway Trolley</b>							
<b>Laryngeal Mask Airway (LMA)</b>							
LMA Size 1 (x1)							
LMA Size 1.5 (x1)							
LMA Size 2 (x1)							
LMA Size 2.5 (x1)							
LMA Size 3 (x1)							
LMA Size 4 (x1)							
<b>Endotracheal Tubes</b>							
ET Tube size 2 uncuffed x1							
ET Tube size 2.5 uncuffed x1							
ET Tube size 3 uncuffed x1							
ET Tube size 3.5 uncuffed x1							
ET Tube size 4 uncuffed x1							
ET Tube size 4.5 uncuffed x1							
ET Tube size 5 uncuffed x1							
ET Tube size 5.5 uncuffed x1							
ET Tube size 6 uncuffed x1							
ET Tube size 6 cuffed x1							
ET Tube size 6.5 cuffed x1							
ET Tube size 7 cuffed x1							
ET Tube size 7.5 cuffed x1							
ET Tube size 8 cuffed x1							
ET Tube size 8.5 cuffed x1							
ET Tube size 8.5 cuffed x1							

Paediatric Resus Checklist

Pink Elastoplast x2 rolls							
Scissors x1							
Durapore x1							
Stylet size 2 (x2)							
Stylet size 4 (x2)							
<b>Airway Accessories</b>							
Catheter mount Adult x1							
Catheter mount Paed x1							
Filter Large (>7.5KG) x1							
Gas Sampling set (Co2) x2							
Magills forceps Adult x1							
Magills forceps Paed x1							
<b>Advanced Airway Kit</b>							
Chryothyroidotomy kit Paediatric x1 (On difficult Airway Trolley Bay 1)							
Tracheostomy tube size 4 (x1)							
Tracheostomy tube size 4.5 (x1)							
Suture pack x1							
4.0 Ethalon reverse cut suture x1							
Mosquito forceps curved x1							
Mosquito forceps straight x1							
Sterile towel x1							
Box of Oxygen bubble tubing x1							
Sterile Jet Insufflation kit x1							
Sterile Scalpel (size 15) x1							
sterile 3way tap with extension x1							
<b>Day</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thur</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>

<b>Insert Date</b>							
Staff Member							
<b>Paediatric Airway Trolley</b>							
<b>Side Basket</b>							
Suction catheters size 6 (x5)							
Suction catheters size 8 (x5)							
Suction catheters size 10 (x5)							
Suction catheters size 12 (x5)							
Tracheal tube boogie Size 5 (x1)							
Tracheal tube boogie Size 10 (x1)							
Tracheal tube boogie Size 15 (x1)							
Stylet extra long x1							
Yankauer suction Paed x5							
Yankauer suction Adult x5							
<b>Paediatric Grab box</b>							
<b>NPA &amp; NG Box (Not sealed)</b>							
NG Tubes Size 8, 10, 12 (x1 Each)							
Bile bag							
Litmus paper							
Catheter tip syringe							
Normal Saline Ampules							
NPA Testing strips							
NPA Trap							

**THE CIRCULATION TROLLEY AND AIRWAY TROLLEY NEED OPENING AND CHECKING ON THE FIRST OF EVERY MONTH – EXPIRY DATES FOR THE NEXT MONTH TAKEN INTO CONSIDERATION IF THIS IS NOT COMPLETED ESCALATE TO NURSE IN CHARGE AND ASK FOR HELP.**