

# Amber Resus Bay 1

Day Insert Date	Mon	Tues	Weds	Thurs	Fri	Sat	Sun
Staff member checking							
Monitoring plugged in and working – SpO2, ECG, BP & Capnography leads							
Disposable kidney dishes							
BP cuffs – Large, Normal, Small adult							
Sharps bin							
Phillips Defib Daily op check, plugged in – ECG check							
Adult defib pads							
ECG electrodes							
Adult bag and valve mask							
Suction unit tubing/filter clean							
X2 suction catheters (Green & White)							
X3 Adult Yankauer							
X1 Bougie 15Fr							
X1 ET stylet 5mm							
Trolley cleaned with tristell							
Blood trolleys cleaned with tristell							
O2 point working None rebreathe mask present							
1000ml and 500ml pressure bags							

# Amber Resus Bay 2 (Paeds)

Day Insert Date	Mon	Tues	Weds	Thurs	Fri	Sat	Sun
Staff member checking							
Monitoring plugged in and working – SpO2, ECG, BP & Capnography leads							
Disposable kidney dishes							
BP cuffs – Large, Normal, Small adult & disposable							
O2 point working with Paed & Adult NRB mask							
Sharps bin							
Phillips Defib Daily op check, plugged in – ECG check							
Adult defib pads							
Paed & neonate defib pads							
ECG electrodes							
Adult bag and valve mask							
1000ml & 500ml pressure bags							
Suction unit tubing/filter clean							
X2 suction catheters (Green & White)							
X3 Adult/Paeds Yankauer							
X1 Bougie 15Fr							
X1 ET stylet 5mm							
Trolley cleaned with tristell							
Blood trolleys cleaned with tristell							

# Amber Resus Daily Checklist

Day Insert Date	Mon	Tues	Weds	Thurs	Fri	Sat	Sun
Staff member checking							
Portable Otoscope (with covers) & Ophthalmoscope **both working**							
Nippy Machine plugged in, stocked & cleaned							
ECG machine: plugged in, stocked and cleaned							
Lucas – check battery is charging and spare suction cups							
Oxylog 3000+ Checked, cleaned and plugged in							
QC BM Machine							
Stock medicines and fluids							
Cannulation trolleys cleaned							
Blood Trolleys cleaned and stocked							
Restock Hygiene and linen cupboard							
Nurses stations clean and tidy							
Donning and doffing room clean & stocked							
Crash trolley checked							

# Amber Resus Grab Boxes

Day Insert Date	Mon	Tues	Weds	Thurs	Fri	Sat	Sun
Staff member checking							
<b>Arterial line</b>							
Arterial Line Box Sealed **SIGN>>>>>>							
Arterial line – Restocked & Sealed (restock list inside box) **SIGN>>>>>>							
<b>Central Line</b>							
Central Line Box Sealed **SIGN>>>>>>							
Central Line Box Restocked & Sealed (restock list inside box) **SIGN>>>>>>							
<b>Difficult Airway</b>							
Difficult Access Box Sealed **SIGN>>>>>>							
Difficult Access Box Restocked & Sealed (restock list inside box) **SIGN>>>>>>							
<b>Chest Drain</b>							
Chest Drain Box Sealed **SIGN>>>>>>							
Chest Drain Box Restocked & Sealed (restock list inside box) **SIGN>>>>>>							
<b>Paediatric Central Line</b>							
Paediatric Central Line Sealed **SIGN>>>>>>							

Paediatric Central Line Restocked & Sealed (restock list inside box) **SIGN>>>>>							
<b>Airway Adjust Box 1</b>							
Airway Adjust Box 1 Sealed **SIGN>>>>>							
Airway Adjust Box 1 Restocked & Sealed (restock list inside box) **SIGN>>>>>							
<b>Airway Adjust Box 2</b>							
Airway Adjust Box 2 Sealed **SIGN>>>>>							
Airway Adjust Box 1 Restocked & Sealed (restock list inside box) **SIGN>>>>>							
<b>Intubation Box 1</b>							
Intubation Box 1 2 Sealed **SIGN>>>>>							
Intubation Box 1 Restocked & Sealed (restock list inside box) **SIGN>>>>>							
<b>Intubation Box 2</b>							
Intubation Box 2 Sealed **SIGN>>>>>							
Intubation Box 2 Restocked & Sealed (restock list inside box) **SIGN>>>>>							
Thoracotomy Tray Sealed & in date							
Amputation Tray Sealed & in date							

SPARE GRAB BOXES KEPT IN MAJAXS ROOM. PLEASE USE TO REFILL THEN REPLACE – THIS WILL HELP WITH THE USING STOCK BEFORE IT EXPIRES.

# Amber Resus Paediatric Trolley

Day Insert Date	Mon	Tue	Wed	Thurs	Fri	Sat	Sun
Staff Member Checking							
<b>Paediatric Airway Trolley</b>							
<b>1<sup>st</sup> Drawer: Oxygen &amp; Basic Airway</b>							
Adult NRB Masks x2							
Paed NRB Mask x2							
2x Child Nebulisers							
Nasal specs neonatal x1							
Nasal specs Paed x1							
Nasal specs Adult x1							
NPA size 3 x1							
NPA size 4 x1							
NPA size 5 x1							
NPA size 6 x1							
NPA size 7 x1							
NPA size 8 x1							
OPA size 000 x2							
OPA size 00 x2							
OPA size 0 x2							
OPA size 1 x2							
OPA size 1.5 x2							
OPA size 2 x2							
OPA size 3 x2							
OPA size 4 x2							
Tongue Depressors x2							
Optilube x2							
Childrens SaO2 probes							
<b>2<sup>nd</sup> Drawer: I-gels and Ventilation</b>							
Face mask size 0							
Face mask size 1							
Face mask size 2							
Face mask size 3							
Face Mask size 4							
Face mask <b>round</b> : size 0							
LMA size 1 x1							
LMA size 1.5 x1							
LMA size 2 x1							
LMA size 2.5 x1							
LMA size 3 x1							
Adult Magill forceps							

Day Insert Date	Mon	Tue	Wed	Thurs	Fri	Sat	Sun
Staff Member Checking							
Paediatric Magills Forceps							
Adult Yankeur (18fr) X2							
Paediatric Yankeur (12fr) x2							
<b>3<sup>rd</sup> Drawer: Intubation</b>							
Laryngoscope handle Adult x1 (Batteries in & working)							
Laryngoscope handle Paed x1 (Batteries in & working)							
Mac 1 x2							
Mac 2 x2							
Mac 3 x1							
Mac 4 x2							
Miller size 00 x1							
Miller size 0 x1							
Miller size 1 x1							
Batteries size C x2							
Batteries size AA x2							
ET Tube cuffed size 6 x1							
ET Tube cuffed size 6.5 x1							
ET Tube cuffed size 7 x1							
ET Tube cuffed size 7.5 x1							
ET Tube cuffed size 8 x1							
ET Tube cuffed size 8.5 x1							
ET Tube uncuffed size 2 x1							
ET Tube uncuffed size 2.5 x1							
ET Tube uncuffed size 3 x1							
ET Tube uncuffed size 3.5 x1							
ET Tube uncuffed size 4 x1							
ET Tube uncuffed size 4.5 x1							
ET Tube uncuffed size 5 x1							
ET Tube uncuffed size 5.5 x1							

Day Insert Date	Mon	Tue	Wed	Thurs	Fri	Sat	Sun
Staff Member Checking							
ET Tube uncuffed size 6 x1							
Adult Catheter Mount							
Paediatric Catheter mount filter							
Capnography lead							
Pink Elastoplast roll							
Sterile Scissors x1							
Durapore x1							
Stylet size 2 x2							
Stylet size 4 x2							
<b>4<sup>th</sup> Drawer: IV Access</b>							
Hypoglycaemia Pack							
IV Packs x3							
Jelco: Yellow x5, Blue x3							
Short extension set x3							
Blue butterfly x3							
Alcohol cleansing wipes x5							
Steristrips x4 small							
Tegaderm x2							
10 x 10 Gauze x3 packs							
5cm K-Band x2							
Micropore roll x1							
5ml Syringe x3							
10ml Syringes x3							
3 way tap x1							
Red filter needles x4							
Blue needles x3							
Blood Culture Bottle (Pink Top) x1							
Blood Bottles FBC, U&E, Glucose, X-Match ***x5 of Each***							
3x Adult VBG syringes							
Blue and Yellow 'Stabbers' x3 each							
Tube of Yellow paraffin							
Capillary Tubes, Caps & Clot Catchers							



Day Insert Date	Mon	Tue	Wed	Thurs	Fri	Sat	Sun
Staff Member Checking							
<b>5<sup>th</sup> Drawer: Intraosseous Access</b>							
EZ-IO Itraosseous Drill							
EZ_IO Needle 3-39 Kg (Pink) x2							
EZ-IO Needle 40 Kg or over (Blue) x2							
EZ-Connect Fluid transfer Set x6 (one set inside each needle box)							
Cook Medical Intraosseous Needle 15.5G x3							
Chloraprep Sponge x2							
Syringe 50ml x1, 10ml x2 & 1ml x2							
Red drawing up needles x2							
Orange needles x2							
0.9% NaCl Flush 10mls x2							
<b>6<sup>th</sup> Drawer: Orogastric Tubes &amp; Water circuits</b>							
OG Tube size 8							
OG Tube size 10							
OG Tube size 12							
Bile Bag							
3 Way Tap x2							
Pen Torch							
10ml Syringe x1							
20ml Syringe x1							
Catheter Tip Syringe x1							
pH Paper							
Adult Mapelson C Circuit Ref 2102000 **unopened and still in packaging**							
Child Mapelson F Circuit *unopened & still in packaging* Ref 2121000							

<b>Day Insert Date</b>	<b>Mon</b>	<b>Tue</b>	<b>Wed</b>	<b>Thurs</b>	<b>Fri</b>	<b>Sat</b>	<b>Sun</b>
Staff Member Checking							
<b>7<sup>th</sup> Drawer: Fluids (Tray not in stacker)</b>							
1x Human Albumen Solution							
2x Giving sets for Alaris SE Pump							
1x Blood giving set for Alaris SE pump							
2x Single Catheter extension sets							
2x 3-Way Taps							
6x Red Drawing up needles							
1x 50ml Leur-Lok Syringe							
6x 20ml Luer-Lok Syringe							
6x 10ml Luer- Lok Syringe							